



SOOO STABLE  
IS **VISALYS<sup>®</sup>** CORE

**CORE BUILD-UP & POST CEMENTATION**



**KETTENBACHDENTAL**

Simply intelligent

# VISALYS® CORE – EASY TREATMENT

## Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



### EXCELLENT STABILITY

- with very good flow properties at the same time



### UNIQUE ACTIVE-CONNECT- TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives



### DENTIN-LIKE GRINDABILITY




### FREE OF BISPHENOL A



Distinguished by the  
experts “The Dental  
Advisor”


# VISALYS® CORE – STABLE RESULTS

## Secure bonding with all relevant adhesives



Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

## Removal depth when grinding various core build-up materials



Visalys® Core can be precisely processed similar to dentin.

# VISALYS® CORE – THE RANGE



## Visalys® Core

<b>Intro pack</b>	White Dentin	1 x 9 g (5 mL) automix syringe, 10 brown mixing tips, 5 intraoral tips, 5 endo tips
<b>Normal pack</b>	White Dentin	2 x 9 g (5 mL) automix syringe, 20 brown mixing tips, 10 intraoral tips, 10 endo tips
<b>Normal pack</b>	White Dentin	1 x 45 g (25 mL) cartridge, 20 yellow mixing tips, 20 intraoral tips



## NOTES

---



---



---



---

Kettenbach SNC  
7, Place de la Gare  
57200 Sarreguemines · France  
info@kettenbach.com  
www.kettenbach-dental.com



022740\_1024

  
**KETTENBACHDENTAL**  
Simply intelligent