



HIGHLY IMPACTFUL – **VISALYS® TEMP**

TEMPORARY CROWN AND BRIDGE MATERIAL



KETTENBACHDENTAL

Simply intelligent

VISALYS® TEMP – ABSOLUTELY BRILLIANT

Temporary crown and bridge material

Visalys® Temp – for the fabrication of temporary crowns, partial crowns, bridges, inlays, onlays and veneers.

EXTREMELY UNBREAKABLE

ALSO SUITABLE AS LONG-TERM
TEMPORARIES

HIGHLY POLISHABLE

EASY & PRECISE PROCESSING


EASY TO REPAIR



Distinguished by the experts
"The Dental Advisor"

VISALYS® TEMP – HIGHLY IMPACTFUL


High breaking strength



Measurements from the University Hospital Regensburg, based on EN ISO 20795-1:2013 (D).

High breaking strength reduces the risk of a temporary denture breaking in the mouth under high or sudden chewing loads

Working and setting times



Note:

Working times apply for products, which have been stored and processed at a temperature of 23°C/74°F and a relative humidity of 50%. Higher temperatures shorten the working time and lower temperatures extend these times.

VISALYS® TEMP – THE RANGE

Visalys® Temp

Normal pack	A1		REF 13780
	A2		REF 13781
	A3	1 x 71 g (50 mL) cartridge,	REF 13782
	A3.5	15 blue-orange mixing tips	REF 13790
	B1		REF 13784
	BL		REF 13788
Bonus pack	A2	5 x 71 g x (50 mL) cartridges,	REF 13794
	A3	15 blue-orange mixing tips	REF 13795



Made
in
Germany

142741_2425



Kettenbach GmbH & Co. KG
Im Heerfeld 7
35713 Eschenburg · Germany
Phone +49 (0)2774 7050
info@kettenbach.com
www.kettenbach-dental.com


KETTENBACHDENTAL
Simply intelligent