



COMFORTABLE UND PRECISE – **IDENTIUM®**

PRECISION IMPRESSION



KETTENBACHDENTAL

Simply intelligent

IDENTIUM[®] – THE EXCEPTIONAL MATERIAL COMBINATION!

A special generation of precision impression material

Combines the advantages of two proven impression materials (A-silicone and polyether) – Vinylsiloxanether[®].




ONE-STEP IMPRESSION

- Double-mix impressions
- Monophase impressions



INDICATIONS

- Crowns and bridges
- Inlays and onlays
- Functional impressions
- Implants
- Combined restorations (fixation impressions)
- Pick-up impressions



Double-mix impression




Implant impression



IDENTIUM® – COMFORTABLE IN ALL DIMENSIONS

Removal from the mouth – as light as a feather

Identium® vs. Impregum™: **2.6 times** lower withdrawal forces



Ease of removal from the patient's mouth avoids distortions and ensures a stress-free impression taking process.


For your precision.



Advantages your patient will appreciate too

Fast final setting time

with sufficient processing time and short intraoral setting time



Sufficient intraoral working time makes it ideal for large cases.
For your precision.

Outstanding flowability



With the outstanding flowability of **Identium® Light**, you can reproduce even the narrowest gaps in great detail. **For your precision.**



- **Odorless and tasteless material** reduces gag reflex
- **Quick setting and easy removal from the mouth** for a relaxed and pleasant experience
- **Initially hydrophilic** – good performance even with residual moisture

COMFORTABLE AND PRECISE IMPRESSIONS



EASY REMOVAL FROM THE MOUTH

SHORT INTRAORAL SETTING TIME

INITIALLY HYDROPHILIC

Available in viscosities
Heavy, Medium and Light



Selectable setting times
regular (4:30 min.)
fast (3:30 min.)



Compatible with most common
mixing devices



IDENTIUM® – THE RANGE

Identium®					
HEAVY	Identium® Heavy		Identium® Heavy Fast		
	Refill pack	REF 14725	REF 14727		
2 × 380 mL cartridges		2 × 380 mL cartridges			
MEDIUM	Identium® Medium		Identium® Medium Fast		
	Refill pack	REF 14717	REF 14719		
2 × 380 mL cartridges		2 × 380 mL cartridges			
LIGHT	Identium® Light		Identium® Light Fast		
	Normal pack	REF 13701	REF 13711		
2 × 50 mL cartridges, 8 grey mixing tips		2 × 50 mL cartridges, 8 grey mixing tips			
ACCESSORIES	Identium® Adhesive	Dynamic Mixers	Mixing tips	Sympress Dispenser	Applyfix 5
	REF 14204	REF 17900	REF 17226 / REF 17240	REF 35910	REF 17204
	10 ml bottle	45 pieces white	50 mixing tips grey / 100 mixing tips grey	Electrical dosing and mixing device	2 application syringes, 12 syringe tips

All 380 mL and 50 mL cartridges are available in a changed cartridge system.



Kettenbach GmbH & Co. KG
 Im Heerfeld 7
 35713 Eschenburg · Germany
 Phone +49 (0)2774 7050
 info@kettenbach.com
 www.kettenbach-dental.com



1410870_0625


KETTENBACHDENTAL
 Simply intelligent